



VISIBLE PATIENT PLANNING

Release version 1.0.17 (April 2022)

This release note represents all changes since the Visible Patient 1.0.15 (Last distributed version).

Software

Major feature

None

New minor features

- **Integrate the following languages:**

- Czech
- Danish
- German
- English
- Spanish
- French
- Italian
- Japanese
- Korean
- Dutch
- Norwegian
- Polish
- Portuguese (Portugal)
- Portuguese (Brazil)
- Slovenian
- Swedish

Fixed bugs

- Fix a software instability when reading DICOM from shared folder.
- Issue with specific characters in DICOM file path.

Remaining problems

None

User manual

- Added Australian sponsor
- Added New Zealand sponsor