



# VISIBLE PATIENT SENDER

*Release version 1.0.17 (April 2022)*

This release note represents all changes since the Visible Patient 1.0.15 (Last distributed version).

## Software

### Major feature

None

### New minor features

- **Integrate the following languages:**

- Czech
- Danish
- German
- English
- Spanish
- French
- Italian
- Japanese
- Korean
- Dutch
- Norwegian
- Polish
- Portuguese (Portugal)
- Portuguese (Brazil)
- Slovenian
- Swedish

## Fixed bugs

- Fix a software instability when reading DICOM from shared folder.
- Issue with specific characters in DICOM file path.

## Remaining problems

None

## User manual

- Added Australian sponsor
- Added New Zealand sponsor