



# VISIBLE PATIENT PLANNING

*Release version 1.0.19 (Dec 2023)*

This release note represents all changes since the Visible Patient 1.0.18 (Last distributed version).

## Software

### Major feature

None

### New minor features

- Improve Organ manager layout on clip applying activity.

### Fixed bugs

- Fix DICOM instability with RGB image.
- Fix DICOM instability with private SOP class instances.
- Fix DICOM instability when reading folder with hidden files.
- Fix software instability with macSonoma OS.

### Remaining problems

None

### User manual

- Update UDI.
- Update link to user manual webpage.
- Update link to authorized countries webpage.

### Labeling

- Update UDI.